

Experiment 6: Infrared Spectroscopy Tutorial

Preparation

There are no prelab questions and the instructor led tutorial during the in-person lab is optional.

Review Sections 12.5-12.8 on interpreting infrared spectra in the course textbook.

The tutorial can be found as an online resource.



Objectives

The purpose of this experiment is to improve your skill in

1. identifying the functional group or groups present in a compound, given a list of the most prominent absorptions in the infrared spectrum and a table of characteristic absorption frequencies.
2. identifying the broad regions of the infrared spectrum to determine the presence of functional groups, such as alcohols, amines, and carbonyl groups, in an unknown compound.

Interpretation of Infrared Spectra

The laboratory instructors will provide a brief infrared tutorial and workshop within the lab class. However, much of this experiment is not wet benchtop chemistry and can be done at home.

1. Review the Theory on Infrared Spectroscopy
2. Review the Listing of Organic Functional Groups and their corresponding Infrared Spectra.
3. Perform the Sample Infrared Spectrum Problems.
4. Answer four (4) of the 'Exp. 6 Infrared Unknown Downloads.'

The Unknown Spectra can be at:



Introduction to Infrared Spectroscopy- Theory and Practice

Electromagnetic Radiation

As you read this page, uncountable numbers of photons or 'light particles' are reflecting off its surface and are being absorbed by pigments (i.e. complex organic molecules) in the rod and cone cells in the retina of your eye. Where the ink (i.e. complex organic dye) has absorbed the photons, you perceive a dark area (i.e. letters) due to the lack of photons from that point on the paper.

On a deeper level, photons (and electrons) are actually wave/particle dualities as described by quantum physics. Photons carry only a discrete amount of energy, called quanta, but the amount of energy of a quanta is defined by the equation, $e = h \nu = h c / l$ where:

- e = the energy of 1 photon (quanta)
- h = Planck's constant (6.62×10^{-27} erg sec)
- ν = Frequency in hertz (cycles or l per sec)
- c = Speed of light (3×10^{10} cm per sec)
- l = Wavelength in cm


Thus, the amount of energy carried by a photon varies directly with its frequency, and because of the relationship between frequency and wavelength, varies inversely with its wavelength. Photons also behave like waves of electromagnetic energy traveling at the speed of light.

Practically speaking however, you need only understand that photons are the messengers that carry the electromagnetic force between electrons and other elementary particles. Electrons, whether free or bound in a covalent bond, are capable of absorbing (or emitting) photons and changing their energy state. This leads to different types of excitation (nuclear transformations, electronic, rotational, nuclear spin changes, bond deformation) depending on the amount of energy carried by the photon. High-energy photons (x-ray, gamma ray, and cosmic ray) can cause ionization of the molecule, while UV photons are involved in electronic interactions. Remember it is the interaction of electrons (via photons) in the outer cloud surrounding atoms that forms the foundation of all chemical reactions.

Infrared Radiation

Infrared radiation is composed of photons with a specific range of wavelengths (7.8×10^{-5} cm to 10^{-2} cm) and frequencies ($\sim 10^{14}$ to 10^{12} Hz). This range includes the near infrared, the infrared and far-infrared regions. The actual wavelengths of interest to most organic chemists are 1.667×10^{-3} cm to 2.5×10^{-4} cm (the 'infrared' region). These wavelengths (λ) are most often expressed as their corresponding wave number (n) where $n = 1/\lambda$, with n measured in cm^{-1} . (e.g. 12.5 to $16.6 \mu\text{m} = 4000$ to 600cm^{-1}).

Infrared carries relatively low levels of energy (e.g. ~ 1 to 10 kcal/mol) which, when absorbed, result in only bond vibrations like stretching and bending, e.g., rocking, scissoring, wagging, and twisting (i.e., bond deformations).








Every molecule, depending on its make-up, is capable of absorbing infrared photons and increasing the intensity of its molecular motions. Different functional groups within the molecule will absorb photons at different infrared wavelengths. Thus, when a spectroscopic wavelength scan is performed on an organic molecule, certain λ will be absorbed while other λ will pass through. Once we have the infrared spectrum of a compound, the spectrum can be analyzed and compared with known infrared absorptions for particular functional groups (see Table 6.1).







The infrared spectrum for a particular molecule can be very complex, consisting of many absorption bands because of the many possible motions each atom can undergo (a non-linear molecule has $3N-6$ normal modes of vibration where N = the number of atoms in the molecule). When analyzing a spectrum, it is important to look at four different regions of the spectrum for the presence or absence of specific absorption peaks. **Note:** you are not required to analyze the fingerprint region.

| Wavenumber cm^{-1} | | | | | |
|-----------------------------|---------------------|---------------------|--|---|--------------------|
| 4000 | 3000 | | 2000 | 1400 | 600 |
| N-H O-H | sp^2 CH | sp^3 CH | $\text{C}\equiv\text{N}$ $\text{C}\equiv\text{C}$ | $\text{C}=\text{C}$ $\text{C}=\text{O}$ $\text{C}=\text{N}$ | fingerprint region |


The following pages contain information to help you understand and interpret infrared spectra.

1. a chart showing the structures of various functional groups, which you need to know.
2. the wavenumbers of the functional groups, to help you locate pertinent absorption bands on an infrared spectrum.
3. Diagrams of the shapes and intensities of various infrared absorption bands, which will help in your interpretation of infrared spectra.
4. Finally, your instructor will lead you through the interpretation of sample infrared spectra representative of various functional groups. Unknown spectra are included to allow you to practice on your own. There is a great deal of information to learn, but the more you practice, the easier it becomes to interpret infrared spectra

| FAMILY NAME | FUNCTIONAL GROUP STRUCTURE | EXAMPLES AND NOMENCLATURE |
|-------------|--|---|
| Alkane |  <p style="text-align: center;">sp³ orbitals</p> | <p>H₃C—CH₃ <i>ethane</i></p>  <i>pentane</i> |
| Alkene |  <p style="text-align: center;">sp² orbitals</p> | <p>H₂C=CH₂ <i>ethene</i></p>  <i>propene</i> |
| Alkyne |  <p style="text-align: center;">sp orbitals</p> | <p>H—C≡C—H <i>ethyne</i> (Acetylene)</p> |

| | | |
|---------|---|---|
| Alcohol |  | $\text{H}_3\text{C}-\text{OH}$ <i>methanol</i>  <i>phenol</i> |
| Ether |  | $\text{H}_3\text{C}-\text{O}-\text{CH}_3$ <i>dimethylether</i>  Tetrahydrofuran |
| Amines | <p>Primary</p>  <p>Secondary</p>  | $\text{H}_3\text{C}-\text{NH}_2$ <i>methylamine</i> $\text{H}_3\text{C}-\text{NH}-\text{CH}_3$ <i>dimethylamine</i> |

Carbonyls:




| | | |
|-----------------|---|---|
| Aldehyde | $\begin{array}{c} \text{O} \\ \\ \text{---C---H} \end{array}$ | $\begin{array}{c} \text{O} \\ \\ \text{H}_3\text{C---C---H} \\ \text{ethanal (Acetaldehyde)} \end{array}$ |
| Ketone | $\begin{array}{c} \text{O} \\ \\ \text{---C---C---C---} \\ \quad \quad \\ \text{---} \quad \text{---} \quad \text{---} \end{array}$ | $\begin{array}{c} \text{O} \\ \\ \text{H}_3\text{C---C---CH}_3 \\ \text{propanone (Acetone)} \end{array}$ |
| Carboxylic Acid | $\begin{array}{c} \text{O} \\ \\ \text{---C---C---OH} \\ \\ \text{---} \end{array}$ | $\begin{array}{c} \text{O} \\ \\ \text{H}_3\text{C---C---OH} \\ \text{ethanoic acid (Acetic acid)} \end{array}$ |
| Ester | $\begin{array}{c} \text{O} \\ \\ \text{---C---C---O---C---} \\ \quad \quad \\ \text{---} \quad \text{---} \quad \text{---} \end{array}$ | $\begin{array}{c} \text{O} \\ \\ \text{H}_3\text{C---C---O---CH}_3 \\ \text{methyl ethanoate (Methyl acetate)} \end{array}$ |
| Amides | $\begin{array}{c} \text{O} \\ \\ \text{---C---N---} \\ \quad \\ \text{---} \quad \text{---} \end{array}$ | $\begin{array}{c} \text{O} \\ \\ \text{H}_3\text{C---C---NH}_2 \\ \text{ethanamide (Acetamide)} \end{array}$ |
| Nitriles | $\begin{array}{c} \text{---} \\ \\ \text{---C---C}\equiv\text{N} \\ \\ \text{---} \end{array}$ | $\begin{array}{c} \text{H}_3\text{C---C}\equiv\text{N} \\ \text{ethanenitrile (Acetonitrile)} \end{array}$ |
| Anhydride | $\begin{array}{c} \text{O} \quad \text{O} \\ \quad \\ \text{---C---O---C---} \\ \diagdown \quad \diagup \\ \text{---} \quad \text{---} \end{array}$ | $\begin{array}{c} \text{O} \quad \text{O} \\ \quad \\ \text{H}_3\text{C---C---O---C---CH}_3 \\ \text{acetic anhydride} \end{array}$ |

Table 6.1 Correlation Table of Infrared Absorption and Functional Group.


| Type of Absorption | Wavenumber (cm ⁻¹) | Intensity of Absorption | Absorption of: |
|-------------------------|--------------------------------|-------------------------|-------------------------------|
| O-H stretch | 3400-3640 | strong, broad | alcohol |
| | 2500-3300 | strong, very broad | carboxylic acid |
| N-H stretch | 3310-3350 | medium ('W' shape) | amine (1°) |
| C-H stretch | 3300 | strong | sp C-H of alkyne |
| | 3030 | medium | aromatic |
| | 3020-3100 | medium | sp ² C-H of alkene |
| | 2850-2960 | medium to strong | sp ³ C-H of alkane |
| | 2750 & 2850 | weak-medium ('W' shape) | O=C-H of aldehyde |
| C≡N stretch | 2210-2260 | medium, sharp | nitrile |
| C≡C stretch | 2100-2260 | medium, sharp | alkyne |
| C=O stretch | 1670-1780 | strong, sharp | carbonyl |
| | 1730-1750 | | ester |
| | 1720-1740 | | aldehyde |
| | 1705-1725 | | ketone |
| | 1700-1725 | | carboxylic acid |
| | 1640-1700 | | amide |
| | ca 1800 and 1760 | | anhydride |
| C=C stretch | 1650-1670 | weak-medium, sharp | alkene |
| | 1600, 1500, 1450 | strong sharp | aromatic |
| C=N stretch | 1640-1670 | medium, sharp | imine |
| N-H bend | 1500-1650 | medium to strong, sharp | amine and amide |
| N=O stretch | 1500-1600 (1540) | strong, sharp | nitro-compound |
| | and 1320-1390 | | |
| C-N stretch | 1030, 1230 | medium | amine |
| C-O stretch | 1050-1150 | strong | alcohol |
| | 1250-1310 | strong broad | ester-conjugated |
| | 1240 | strong, broad | ester-acetates |
| | 1175 | strong, broad | ester-unconjugated |
| C-Cl stretch (terminal) | 600-800 | strong | alkyl halide |
| Ar-Cl stretch | 1000-1175 | medium-strong | aryl halide |
| C-Br stretch (terminal) | 500-760 | strong | alkyl halide |
| C-I (terminal) | 500 | strong | alkyl halide |

Note: when a C=C bond is in conjugation with a carbonyl, the observed carbonyl absorption frequency will be $< \sim 30$ cm⁻¹.

Shapes of Infrared Absorption Bands Observed for Different Functional Groups



Typical Absorption Band Shapes (cont.)



How to Interpret an Infrared Spectrum

Step 1 Divide the infrared spectrum into four main areas (use pencil and ruler and take into account any off-shift in the spectrum's wavenumbers).

- i) Above 3000 cm^{-1}
- ii) Between 3000 and 2000 cm^{-1}
- iii) Between 2000 and 1400 cm^{-1}
- iv) Below 1400 cm^{-1} (fingerprint region)

Step 2 Starting at the left of the spectrum, examine the area **above 3000 cm^{-1}** , first looking in the region near 3300 cm^{-1} and record in tabular format the presence/absence of:

- i) a broad, very strong absorption band of an '**O-H**'. If present, it means you know that your molecule is at least an **alcohol**.
- ii) A broad, weak to medium strength, double or single absorption band of '**N-H**'. If present it means you have an **amine** (1° or 2°) or possibly an **amide**.
- iii) A sharp, medium to strong, single absorption band of ' **$\equiv\text{C-H}$** ' of a **terminal alkyne**.

Note: If present, it means you should also see a ' **$\text{C}\equiv\text{C}$** ' absorption near 2250 cm^{-1} .

After examining the region around 3300 cm^{-1} , look for any sharp, weak to medium absorption just above 3000 cm^{-1} (e.g., 3050 cm^{-1}) resulting from the '**C-H**' stretch of a sp^2 hybridized carbon. If present, it means you have a ' **$\text{C}=\text{C-H}$** ' of an alkene or aromatic compound.

Step 3 Next examine the area between **3000 and 2000 cm^{-1}** and record the presence/absence of absorption bands or peaks.

- i) First look just below 3000 cm^{-1} (e.g., $2850\text{-}2950\text{ cm}^{-1}$) resulting from the '**C-H**' stretch of a sp^3 hybridized carbon. If present, it means you are seeing the '**C-H**' stretch of an **$-\text{CH}_2$ or $-\text{CH}_3$** group. Note: This absorption is not very informative as most organic compounds have $-\text{CH}_2$ or $-\text{CH}_3$ groups.
- ii) Then look for the extremely broad peak, actually starting at 3300 cm^{-1} and extending all the way to $\sim 2500\text{ cm}^{-1}$, caused by the **O-H dimer** between two **carboxylic acid** molecules (COOH). This absorption is probably the most difficult to see as other absorption peaks may be overlapping the broad peak.
- iii) Finally look for a sharp, weak to medium peak caused by either ' **$\text{C}\equiv\text{C}$** ' or ' **$\text{C}\equiv\text{N}$** '.
- iv) If present, then the compound is an alkyne (might also have the '**C-H**' of a terminal alkyne, see step 2 above) or a nitrile.

Step 4 Next examine the area between **2000 and 1400 cm^{-1}** and record the presence/absence of absorption bands or peaks.

- i) First look near 1700 cm^{-1} (e.g. $1680\text{-}1750\text{ cm}^{-1}$) for a sharp, strong peak resulting from the '**C=O**' stretch of a **carbonyl**. Note: This absorption is very informative and will be present if your compound is an aldehyde, ketone, ester, amide, or carboxylic acid.
- ii) Next look near 1650 cm^{-1} (e.g. $1600\text{-}1670\text{ cm}^{-1}$) for a sharp, weak peak resulting from the '**C=C**' stretch of an **alkene**.
- iii) Finally look near 1600 cm^{-1} and 1500 cm^{-1} for a sharp, double peak resulting from the '**C=C**' stretch of an **aromatic ring**.

Step 5 If you dare, you may look in the **fingerprint region (area below 1400 cm^{-1})** and record the presence of absorption bands or peaks.

- i) First look near 1200 ($1160\text{-}1310$) cm^{-1} for a sharp, strong peak resulting from the '**C-O**' stretch of an **ester**.
Note: This absorption is very difficult to see and may or may not be present, i.e. conclusive if present, inconclusive if not present.
- ii) If you suspect you have an aromatic ring (absorption bands at ~ 3030 and 1600 and 1500 cm^{-1} present), you may try to discern the substitution pattern of the benzene ring by looking at the strong absorption bands of the **ring 'C-H'** out-of-plane bending vibrations in the region $680\text{-}900\text{ cm}^{-1}$.

| Benzene Substitution Pattern | Ring 'C-H' Absorption Bands Present (cm^{-1}) |
|------------------------------|--|
| monosubstituted | 2 sharp peaks, 730-770, 690-710 |
| <i>ortho</i> disubstituted | 1 sharp peak, 735-770 |
| <i>meta</i> disubstituted | 3 sharp peaks, 860-900, 750-810, 680-725 |
| <i>para</i> disubstituted | 1 sharp peak, 800-860 |
| 1,2,3 trisubstituted | 2 sharp peaks, 760-780, 705-745 |
| 1,3,5 trisubstituted | 2 sharp peaks, 810-865, 675-730 |
| 1,2,4 trisubstituted | 2 sharp peaks, 870-885, 805-825 |

Ref: McMurry, J., 1992. Organic Chemistry, 3rd ed, Brooks/Cole, p.549-550, (4th ed, p.559)
Nakanishi, K., 1964. Infrared Absorption Spectroscopy, Holden Day p.27.

- iii) Again, if you have an aromatic, you may also try to discern the ring substitution pattern of the benzene ring by looking at the very weak overtone-combination absorption bands of the **ring 'C-H'** stretch vibrations in the region $1670\text{-}2000\text{ cm}^{-1}$.

| Benzene Substitution Pattern | Ring 'C-H' Overtone Bands Present (cm^{-1}) |
|------------------------------|--|
| monosubstituted | 4 weak equally spaced and shaped sharp peaks |
| <i>ortho</i> disubstituted | 3 weak irregularly spaced/shaped sharp peaks |
| <i>meta</i> disubstituted | 2 weak sharp peaks + one weak broad peak |
| <i>para</i> disubstituted | 2 weak sharp peaks |

- iv) If you suspect you have a long straight chain (>4 C) alkane, (absorption bands at $2850\text{-}2950\text{ cm}^{-1}$ present but not much else), you may try to see the sharp, weak absorption due to the concerted rocking of >4 -CH_2 in a chain. It lies in the region $720 \pm 10\text{ cm}^{-1}$.

Step 6 Finally, you will summarize your results by making a statement about what functional groups you suspect to be present in the molecule or perhaps you will be asked to select from a list of suggested structures, which molecule most likely would generate the spectrum just analyzed.

Instructor Led Group Infrared Analysis Problems

Use the tables below to record your results of the 'Infrared Spectral Analyses' for the following compounds (infrared spectra appear on the following pages of this lab manual). Label the absorption bands.

| Cyclohexanol | Absorption Band Code# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|----------------------------|-----------------------|--------------------------------|---------------------------|---|-----------------------------|
| >3000 cm ⁻¹ | 1 | 3331 | broad | strong | O-H stretch alcohol |
| 3000-2000 cm ⁻¹ | 2 | 2932 & 2855 | sharp | strong | C-H sp ³ stretch |
| 2000-1500 cm ⁻¹ | none | | | | |
| (Fingerprint) | 3 | 1068 | sharp | strong | C-O of alcohol |

Functional Group absent: no ≡C-H, no N-H, no sp² H-C=, no C≡C, no C≡N, no C=O, no C=C alkene or aromatic

| 2-methyl-3-butyn-2-ol | Absorption Band Code# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|------------------------------|-----------------------|--------------------------------|---------------------------|---|----------------------------|
| >3000 cm ⁻¹ | 1 | ~3380 | broad | strong | O-H stretch alcohol |
| | 2 | 3303 | sharp | strong | |
| 3000-2000 cm ⁻¹ | 3 | 2876,2938,2987 | sharp | med-str. | |
| | 4 | 2120 | sharp | weak | |
| 2000-1500 cm ⁻¹ | none | | | | |

Functional Group absent: no N-H, no sp² H-C=, no C≡N, no C=O, no C=C alkene or aromatic

| 3-buten-2-ol | Absorption Band Code# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|----------------------------|-----------------------|--------------------------------|---------------------------|---|----------------------------|
| >3000 cm ⁻¹ | 1 | ~3350 | broad | | |
| | 2 | 3083 & 3012 | | strong | C-H stretch |
| 3000-2000 cm ⁻¹ | 3 | | sharp | | C-H stretch |
| 2000-1500 cm ⁻¹ | 4 | 1646 | | | |

Functional Group absent: no ≡C-H, no N-H, no C≡C, no C≡N, no C=O, no C=C aromatic

| benzhydrol | Absorption Band Code# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|----------------------------|-----------------------|--------------------------------|---------------------------|---|----------------------------|
| >3000 cm ⁻¹ | 1 | 3392-3359 | broad | | |
| | 2 | 3049 & 3027 | sharp | | C-H stretch |
| 3000-2000 cm ⁻¹ | 3 | 2900 | sharp | | C-H stretch |
| 2000-1500 cm ⁻¹ | 4 | 1598,1495,1458 | sharp | | |

Functional Group absent: no ≡C-H, no N-H, no C≡C, no C≡N, no C=O, no C=C alkene

Instructor Led Group Infrared Analysis Problems (cont.)

| benzaldehyde | Absorption Band Code# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|---------------------|-----------------------|--------------------------------|---------------------------|---|----------------------------|
| | | | | | |
| | | | | | |
| | | | | | |

Functional Group absent: no O-H, no ≡C-H, no N-H, no sp³ C-H, no C≡C, no C≡N, no C=C alkene

| acetic acid | Absorption Band Code# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|--------------------|-----------------------|--------------------------------|---------------------------|---|----------------------------|
| | | | | | |
| | | | | | |
| | | | | | |


Functional Group absent:

| dibutylamine | Absorption Band Code# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|---------------------|-----------------------|--------------------------------|---------------------------|---|----------------------------|
| | | | | | |
| | | | | | |
| | | | | | |


Functional Group absent:

CHEM350 Lab Manual 2024

| | | | |
|--------------|------------|-------------|-----------------------------|
| HIT-NO=1077 | SCORE= () | SDBS-NO=581 | IR-NIDA-09018 : LIQUID FILM |
| CYCLOHEXANOL | | | |
| $C_6H_{12}O$ | | | |




| | | | | | | | | | |
|------|----|------|----|------|----|-----|----|-----|----|
| 3331 | 13 | 1704 | 86 | 1256 | 47 | 970 | 21 | 667 | 64 |
| 2932 | 4 | 1467 | 42 | 1238 | 53 | 926 | 84 | 557 | 53 |
| 2855 | 8 | 1452 | 17 | 1174 | 86 | 890 | 39 | 462 | 61 |
| 2696 | 68 | 1363 | 29 | 1140 | 55 | 863 | 81 | | |
| 2666 | 62 | 1346 | 41 | 1068 | 11 | 845 | 47 | | |
| 2568 | 74 | 1329 | 50 | 1034 | 52 | 835 | 74 | | |
| 2233 | 84 | 1298 | 48 | 1025 | 32 | 789 | 64 | | |




CHEM350 Lab Manual 2024

| | | | |
|---------------------------------|------------|--------------|-----------------------------|
| HIT-NO=2057 | SCORE= () | SDBS-NO=2818 | IR-NIDA-00603 : LIQUID FILM |
| 2-METHYL-3-BUTYN-2-OL | | | |
| C ₅ H ₈ O | | | |




| | | | | | |
|------|----|------|----|------|----|
| 3939 | 81 | 2120 | 84 | 1170 | 8 |
| 3303 | 4 | 1706 | 74 | 963 | 10 |
| 2987 | 7 | 1465 | 52 | 889 | 22 |
| 2938 | 38 | 1459 | 49 | 708 | 42 |
| 2876 | 62 | 1378 | 21 | 652 | 20 |
| 2811 | 77 | 1368 | 17 | 558 | 30 |
| 2619 | 81 | 1216 | 18 | | |




CHEM350 Lab Manual 2024


| | | | |
|---------------------------------|------------|--------------|-----------------------------|
| HIT-NO=1306 | SCORE= () | SDBS-NO=1211 | IR-NIDA-01713 : LIQUID FILM |
| 3-BUTEN-2-OL | | | |
| C ₄ H ₈ O | | | |



| | | | | | | | |
|------|----|------|----|------|----|-----|----|
| 3083 | 33 | 1646 | 62 | 1066 | 12 | 686 | 41 |
| 3012 | 36 | 1452 | 33 | 1005 | 24 | 539 | 70 |
| 2979 | 4 | 1424 | 20 | 991 | 15 | 443 | 84 |
| 2932 | 26 | 1371 | 29 | 978 | 38 | | |
| 2876 | 24 | 1292 | 35 | 922 | 5 | | |
| 1846 | 79 | 1152 | 29 | 829 | 80 | | |
| 1711 | 81 | 1103 | 14 | 810 | 72 | | |



| | | | |
|-----------------|------------|-------------|--------------------------|
| HIT-NO=1144 | SCORE= () | SDBS-NO=869 | IR-NIDA-47684 : KBR DISC |
| BENZHYDROL | | | |
| $C_{13}H_{12}O$ | | | |




| | | | | | | | | | |
|------|----|------|----|------|----|------|----|-----|----|
| 3392 | 26 | 1698 | 72 | 1346 | 66 | 1172 | 44 | 863 | 72 |
| 3359 | 30 | 1495 | 17 | 1317 | 64 | 1156 | 72 | 754 | 16 |
| 3331 | 29 | 1458 | 21 | 1289 | 86 | 1037 | 28 | 735 | 15 |
| 3086 | 64 | 1446 | 27 | 1271 | 39 | 1032 | 22 | 702 | 4 |
| 3059 | 57 | 1422 | 70 | 1244 | 72 | 1019 | 11 | 654 | 37 |
| 3049 | 70 | 1395 | 72 | 1201 | 88 | 1003 | 60 | 604 | 39 |
| 3027 | 63 | 1361 | 48 | 1182 | 46 | 912 | 72 | 646 | 62 |


Oc1ccccc1-c2ccccc2

CHEM350 Lab Manual 2024


| | | | |
|---------------------------------|------------|-------------|-------------------------------|
| HIT-NO=1100 | SCORE= () | SDBS-NO=672 | IR-NIDA-08667 : CCL4 SOLUTION |
| BENZALDEHYDE | | | |
| C ₇ H ₆ O | | | |




| | | | | | | | |
|------|----|------|----|------|----|------|----|
| 3088 | 86 | 1778 | 86 | 1388 | 81 | 1024 | 84 |
| 3068 | 86 | 1709 | 4 | 1310 | 88 | 828 | 55 |
| 3032 | 86 | 1653 | 84 | 1285 | 86 | 688 | 47 |
| 2847 | 86 | 1607 | 81 | 1203 | 39 | 673 | 86 |
| 2814 | 79 | 1598 | 72 | 1167 | 88 | 650 | 64 |
| 2732 | 77 | 1588 | 72 | 1160 | 84 | | |
| 2693 | 86 | 1466 | 74 | 1072 | 86 | | |



| | | | |
|--------------|------------|--------------|-----------------------------|
| HIT-NO=1506 | SCORE= () | SDBS-NO=1672 | IR-NIDA-04481 : LIQUID FILM |
| DIBUTYLAMINE | | | |
| $C_8H_{19}N$ | | | |



| | | | | | |
|------|----|------|----|------|----|
| 3684 | 86 | 1466 | 20 | 1181 | 79 |
| 3290 | 72 | 1378 | 43 | 1131 | 26 |
| 2959 | 4 | 1338 | 70 | 1096 | 84 |
| 2929 | 5 | 1302 | 72 | 966 | 79 |
| 2874 | 9 | 1245 | 81 | 961 | 81 |
| 2862 | 10 | 1234 | 79 | 898 | 79 |
| 2810 | 16 | 1206 | 79 | 736 | 41 |
| 1617 | 86 | 892 | 41 | | |



Infrared Analysis Practice Problems

Use the tables below to record your results of the 'Infrared Spectral Analyses' of the provided known spectra in this lab manual.

| benzaldehyde | Absorption Band# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|---------------------|------------------|--------------------------------|---------------------------|---|----------------------------|
| | | | | | |
| | | | | | |
| | | | | | |

Functional Group(s) absent:

| benzoic acid | Absorption Band# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|---------------------|------------------|--------------------------------|---------------------------|---|----------------------------|
| | | | | | |
| | | | | | |
| | | | | | |

Functional Group(s) absent:


| phenylacetylene | Absorption Band# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|------------------------|------------------|--------------------------------|---------------------------|---|----------------------------|
| | | | | | |
| | | | | | |
| | | | | | |

Functional Group(s) absent:


| styrene | Absorption Band# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|----------------|------------------|--------------------------------|---------------------------|---|----------------------------|
| | | | | | |
| | | | | | |
| | | | | | |

Functional Group(s) absent:

| | | | |
|---------------------------------|------------|-------------|-------------------------------|
| HIT-NO=1100 | SCORE= () | SDBS-NO=672 | IR-NIDA-08667 : CCL4 SOLUTION |
| BENZALDEHYDE | | | |
| C ₇ H ₆ O | | | |




| | | | | | | | |
|------|----|------|----|------|----|------|----|
| 3088 | 86 | 1778 | 86 | 1388 | 81 | 1024 | 84 |
| 3068 | 86 | 1709 | 4 | 1310 | 88 | 828 | 55 |
| 3032 | 86 | 1653 | 84 | 1285 | 86 | 688 | 47 |
| 2847 | 86 | 1607 | 81 | 1203 | 39 | 673 | 86 |
| 2814 | 79 | 1598 | 72 | 1167 | 88 | 650 | 64 |
| 2732 | 77 | 1588 | 72 | 1160 | 84 | | |
| 2693 | 86 | 1466 | 74 | 1072 | 86 | | |
| 2654 | 84 | 1634 | 68 | 1266 | 64 | 909 | 36 |




HIT-NO=1451 SCORE= () SDBS-NO=1460 IR-NIDA-04316 : LIQUID FILM

ETHYL BENZOATE


$C_9H_{10}O_2$




| | | | | | | | | | |
|------|----|------|----|------|----|------|----|-----|----|
| 3091 | 81 | 1959 | 84 | 1466 | 58 | 1109 | 12 | 607 | 61 |
| 3064 | 72 | 1822 | 84 | 1452 | 23 | 1071 | 26 | 782 | 79 |
| 3035 | 77 | 1719 | 6 | 1392 | 57 | 1029 | 21 | 711 | 9 |
| 2983 | 46 | 1603 | 49 | 1368 | 31 | 1002 | 60 | 688 | 47 |
| 2939 | 88 | 1585 | 53 | 1315 | 23 | 937 | 79 | 675 | 52 |
| 2907 | 70 | 1492 | 72 | 1276 | 4 | 873 | 74 | 618 | 61 |
| 2874 | 79 | 1478 | 64 | 1176 | 34 | 861 | 70 | 606 | 84 |



| | | | |
|--|------------|-------------|--------------------------|
| HIT-NO=1081 | SCORE= () | SDBS-NO=673 | IR-NIDA-63340 : KBR DISC |
| BENZOIC ACID | | | |
| C ₇ H ₆ O ₂ | | | |




| | | | | | | | | | |
|------|----|------|----|------|----|------|----|-----|----|
| 3073 | 49 | 2678 | 57 | 1426 | 35 | 1112 | 77 | 936 | 42 |
| 3012 | 53 | 2607 | 62 | 1327 | 17 | 1107 | 77 | 812 | 62 |
| 2996 | 53 | 2564 | 57 | 1294 | 14 | 1102 | 74 | 805 | 70 |
| 2986 | 53 | 1689 | 5 | 1187 | 58 | 1074 | 62 | 708 | 4 |
| 2886 | 52 | 1603 | 58 | 1180 | 60 | 1028 | 53 | 685 | 55 |
| 2836 | 52 | 1585 | 59 | 1129 | 84 | 1001 | 72 | 667 | 52 |
| 2726 | 66 | 1464 | 28 | 1118 | 77 | 943 | 50 | 654 | 68 |




HIT-NO=1445 SCORE= () SDBS-NO=1444 IR-NIDA-63379 : LIQUID FILM

PHENYLACETYLENE


C₈H₆



| | | | | | | | | | |
|------|----|------|----|------|----|------|----|-----|----|
| 3906 | 86 | 2110 | 74 | 1674 | 68 | 1176 | 81 | 843 | 86 |
| 3306 | 16 | 1954 | 84 | 1488 | 15 | 1159 | 81 | 757 | 4 |
| 3291 | 12 | 1900 | 86 | 1444 | 31 | 1100 | 84 | 692 | 9 |
| 3081 | 62 | 1886 | 84 | 1386 | 86 | 1071 | 67 | 666 | 24 |
| 3058 | 62 | 1808 | 86 | 1331 | 81 | 1026 | 53 | 621 | 23 |
| 3034 | 66 | 1757 | 86 | 1262 | 72 | 918 | 57 | 530 | 30 |
| 3022 | 72 | 1698 | 70 | 1246 | 66 | 883 | 84 | 514 | 43 |




| | | | |
|---------------------------------|------------|-------------|-----------------------------|
| HIT-NO=1114 | SCORE= () | SDBS-NO=669 | IR-NIDA-05064 : LIQUID FILM |
| BENZONITRILE | | | |
| C ₇ H ₅ N | | | |



| | | | | | | | | | |
|------|----|------|----|------|----|------|----|------|----|
| 3646 | 86 | 2178 | 86 | 1682 | 72 | 1288 | 60 | 1001 | 74 |
| 3412 | 86 | 1969 | 84 | 1572 | 81 | 1193 | 70 | 927 | 60 |
| 3088 | 74 | 1899 | 81 | 1492 | 17 | 1178 | 66 | 846 | 64 |
| 3066 | 62 | 1816 | 84 | 1448 | 12 | 1163 | 77 | 768 | 4 |
| 3004 | 84 | 1768 | 86 | 1441 | 72 | 1098 | 81 | 688 | 6 |
| 2256 | 84 | 1688 | 84 | 1392 | 84 | 1072 | 64 | 625 | 64 |
| 2230 | 8 | 1699 | 66 | 1336 | 81 | 1027 | 66 | 648 | 8 |

| | | | |
|-------------------------------|------------|--------------|-----------------------------|
| HIT-NO=2170 | SCORE= () | SDBS-NO=3044 | IR-NIDA-10290 : LIQUID FILM |
| STYRENE | | | |
| C ₈ H ₈ | | | |



| | | | | | | | | | |
|------|----|------|----|------|----|------|----|-----|----|
| 3299 | 74 | 1946 | 74 | 1496 | 9 | 1202 | 68 | 841 | 74 |
| 3082 | 28 | 1876 | 74 | 1449 | 23 | 1182 | 74 | 777 | 5 |
| 3060 | 29 | 1821 | 66 | 1412 | 43 | 1156 | 77 | 738 | 72 |
| 3027 | 27 | 1689 | 77 | 1383 | 77 | 1083 | 43 | 698 | 4 |
| 3009 | 38 | 1630 | 24 | 1334 | 70 | 1021 | 33 | 650 | 70 |
| 2980 | 57 | 1601 | 59 | 1317 | 86 | 992 | 14 | 555 | 60 |
| 1965 | 77 | 1576 | 41 | 1290 | 66 | 909 | 9 | 442 | 67 |

Write-up

Fill in the following form below and answer the post-lab questions. Use the Word version of the report form so you can add additional space for your answers. You will also need to download four (4) unknown spectra and include that in your report. When complete save as a PDF and email as an attachment to your Academic Expert for grading.

CHEM 350 Experiment 6 Report Form

Infrared Spectroscopy Tutorial

Date: _____

Student Name: _____ ID Number: _____

Infrared Knowns

Fill in the following three (3) analyses tables to reflect your characterization of the spectra provided (above).

| | Absorption Band# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|----------------------|------------------|--------------------------------|---------------------------|---|----------------------------|
| cyclohexanone | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Functional Group(s) absent:

| | Absorption Band# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|-----------------------|------------------|--------------------------------|---------------------------|---|----------------------------|
| ethyl benzoate | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Functional Group(s) absent:

| | Absorption Band# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|-------------------|------------------|--------------------------------|---------------------------|---|----------------------------|
| benzotrile | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Functional Group(s) absent:

Infrared Unknowns

Select four (4) unknowns from the 'Exp. 6 Infrared Unknown Downloads' list:



Download 4 of the possible 20 spectra (PDFs). Please neatly fill out the table on the unknown spectra and remember to fully label each of the absorption bands identified and identify the compound. If you find the tables on the PDFs too small use this Word template to give yourself more space to write/type.

| Code: Name: | Absorption Band# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|----------------|---------------------|-----------------------------------|------------------------------------|--|----------------------------------|
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Functional Group absent:

| Code: Name: | Absorption Band# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|----------------|---------------------|-----------------------------------|------------------------------------|--|----------------------------------|
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Functional Group absent:

| Code: Name: | Absorption Band# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|----------------|---------------------|-----------------------------------|------------------------------------|--|----------------------------------|
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |



Functional Group absent:

| Code: Name: | Absorption Band# | Wavenumber (cm ⁻¹) | Peak Shape (sharp, broad) | Peak Intensity (strong, medium or weak) | Functional Group Indicated |
|----------------|---------------------|-----------------------------------|------------------------------------|--|----------------------------------|
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Functional Group absent:

Questions

1. What are the major differences you would see in the infrared spectra of an alkane, alkene, and alkyne?
2. Consider the C=O absorption of three compounds: 2-butanone (1715 cm^{-1}), propanoyl chloride (1772 cm^{-1}), and propyl amide (1650 cm^{-1}). Explain the observed differences.
3. Describe how IR spectroscopy might be used to monitor the progress of each of the following reactions.



I certify that this submitted laboratory report represents entirely my own efforts. I have read and understand the Athabasca University policies regarding, and sanctions for, plagiarism.

Signature: _____ Date: _____